

Title	Semi-Volatile Organic Compounds
Specimen Requirements	<p>Samples for SVOC analysis will only be accepted from the Indiana State Department of Health and other government agencies.</p> <p>Contact: Mark Starzynski 921-5580 for sample containers and instructions.</p>
Sampling Materials	Water must be collected in one liter glass bottles with teflon lined caps, both of which were previously rinsed with solvents. Solids are collected in solvent rinsed pint glass jars with clean aluminum foil between the jar and the cap.
Procedural Notes	<p>Water --- fill 2 of the one liter glass bottles to the top.</p> <p>Solid ---- fill 1 or 2 of the pint glass jars with the solid.</p>
Shipping Instructions	<p>Samples must be kept refrigerated to 4 °C from the time of collection to the time of analysis. All samples must be extracted within seven days of collection and completely analyzed within forty days of extraction.</p> <p>Indiana State Department of Health Chemistry Laboratories 550 W. 16th Street, Suite B Indianapolis, In 46202</p>
Reporting and TAT	Results will be e-mailed to the submitter normally within 30 days.